

Product Information	
Catalog No.	RA14535
Product name	Anti-Fel d 1 Human IgE Antibody (11A12)
Applications	ELISA
Species reactivity	Felis catus (Cat) (Felis silvestris catus)
Host species	Human
Clonality	Monoclonal
Isotype	IgE
Clone ID	11A12
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Target	Fel d 1, Felis domesticus allergen.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.