

Product Information	
Catalog No.	RA14553
Product name	Anti-Ara h 3 Human IgE Antibody (3C3)
Applications	ELISA
Species reactivity	Arachis hypogaea (Peanut)
Host species	Human
Clonality	Monoclonal
Isotype	IgE
Clone ID	3C3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Target	Ara h 3, Arachis hypogaea allergen.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.