


Product Information	
Catalog No.	RA14564
Product name	InVivoMAb Anti-Human CD9 Antibody (Iv0274)
Applications	ELISA, FuncS
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0274
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Motility-related protein, Leukocyte antigen MIC3, TSPAN29, Cell growth-inhibiting gene 2 protein, CD9, MIC3, p24, Tspan-29, MRP-1, 5H9 antigen, Tetraspanin-29, CD9 antigen
Accession	P21926
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human CD9 Antibody (lv0274)

SDS-PAGE