


Product Information	
Catalog No.	RA14565
Product name	InVivoMAb Anti-HMPV F/Fusion glycoprotein F0 Antibody (MPV467)
Applications	ELISA, Neutralization
Species reactivity	Human metapneumovirus (strain CAN97-83) (HMPV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	MPV467
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Fusion glycoprotein F0, Protein F, Fusion glycoprotein F2, Fusion glycoprotein F1
Accession	Q6WB98
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-HMPV F/Fusion glycoprotein F0 Antibody (MPV467)

SDS-PAGE