


Product Information	
Catalog No.	RA14578
Product name	Anti-Human CXCL4/PF4 Antibody (1E12#)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	1E12#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	C-X-C motif chemokine 4, PF-4, Oncostatin-A, Endothelial cell growth inhibitor, CXCL4, Iroplact, Platelet factor 4, PF4, SCYB4
Accession	P02776
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CXCL4/PF4 Antibody (1E12#)

SDS-PAGE