


Product Information	
Catalog No.	RA14589
Product name	Anti-Human RIOX1/C14ORF169 IP Antibody (SAA0194)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0194
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	ROX, MAPJD, C14orf169, hsNO66, Myc-associated protein with JmjC domain, Ribosomal oxygenase 1, Nucleolar protein 66, Ribosomal oxygenase NO66, Histone lysine demethylase NO66, 60S ribosomal protein L8 histidine hydroxylase, NO66, Bifunctional lysine-specific demethylase and histidyl-hydroxylase NO66, RIOX1
Accession	Q9H6W3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human RIOX1/C14ORF169 IP Antibody (SAA0194)

SDS-PAGE