


Product Information	
Catalog No.	RA14596
Product name	Anti-Human KDM2A IP Antibody (SAA0201)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0201
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	KDM2A, CXXC-type zinc finger protein 8, FBXL11, JHDM1A, FBL7, KIAA1004, F-box protein FBL7, Lysine-specific demethylase 2A, FBL11, [Histone-H3]-lysine-36 demethylase 1A, F-box protein Lilina, JmjC domain-containing histone demethylation protein 1A, F-box/LRR-repeat protein 11, CXXC8, F-box and leucine-rich repeat protein 11
Accession	Q9Y2K7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human KDM2A IP Antibody (SAA0201)

SDS-PAGE