


Product Information	
Catalog No.	RA14604
Product name	Anti-Human L3MBTL3 IP Antibody (SAA0209)
Applications	ELISA, IF, IP, IP-MS
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0209
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Lethal(3)malignant brain tumor-like protein 3, H-l(3)mbt-like protein 3, KIAA1798, L3MBTL3, L3mbt-like 3, MBT-1, L(3)mbt-like protein 3, MBT1
Accession	Q96JM7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human L3MBTL3 IP Antibody (SAA0209)

SDS-PAGE