


Product Information	
Catalog No.	RA14611
Product name	Anti-Human PHF8 IP Antibody (SAA0216)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0216
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	[histone H3]-dimethyl-L-lysine(9) demethylase PHF8, KIAA1111, [histone H3]-dimethyl-L-lysine(36) demethylase PHF8, PHF8, ZNF422, Histone lysine demethylase PHF8, PHD finger protein 8
Accession	Q9UPP1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human PHF8 IP Antibody (SAA0216)

SDS-PAGE