


Product Information	
Catalog No.	RA14621
Product name	Anti-Human SFMBT2 IP Antibody (SAA0226)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0226
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SFMBT2, KIAA1617, Scm-like with 4 MBT domains protein 2, Scm-like with four MBT domains protein 2
Accession	Q5VUG0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human SFMBT2 IP Antibody (SAA0226)

SDS-PAGE