


Product Information	
Catalog No.	RA14643
Product name	Anti-Human RPRD1A IP Antibody (SAA0248)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0248
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Regulation of nuclear pre-mRNA domain-containing protein 1A, p15INK4B-related protein, Cyclin-dependent kinase inhibitor 2B-related protein, RPRD1A, P15RS
Accession	Q96P16
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human RPRD1A IP Antibody (SAA0248)

SDS-PAGE