


Product Information	
Catalog No.	RA14646
Product name	Anti-Human DNMT1 IP Antibody (SAA0251)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0251
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CXXC-type zinc finger protein 9, MCMT, AIM, DNMT1, DNA MTase HsaI, DNMT, DNA (cytosine-5)-methyltransferase 1, Dnmt1, DNA methyltransferase HsaI, M.HsaI, CXXC9
Accession	P26358
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human DNMT1 IP Antibody (SAA0251)

SDS-PAGE