


Product Information	
Catalog No.	SC0006
Product name	Anti-Mouse IgG1 Fc Antibody (SAA0489)
Applications	ELISA, IHC, IF, IP, WB
Species reactivity	Mouse
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	SAA0489
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ig gamma-1 chain C region secreted form, Ighg1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse IgG1 Fc Antibody (SAA0489)

SDS-PAGE