


Product Information	
Catalog No.	SC0018
Product name	Anti-Biotin Antibody (SAA0349)
Applications	SPR, WB
Species reactivity	General
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA0349
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Biotin, Vitamin B7, Vitamin H, BTN
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Biotin Antibody (SAA0349)

SDS-PAGE