

Product Information	
Catalog No.	SC0027
Product name	Anti-Human IgM Monoclonal Antibody (1A185)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	1A185
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	Second antibody, Second Ab, Mouse Anti-Human IgM Monoclonal Antibody
Immunogen	IgM
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.