

| Product Information   |                                                                |
|-----------------------|----------------------------------------------------------------|
| Catalog No.           | SC0033                                                         |
| Product name          | Anti-Human IgM Monoclonal Antibody (1A184), HRP                |
| Applications          | ELISA, WB                                                      |
| Species reactivity    | Human                                                          |
| Host species          | Mouse                                                          |
| Clonality             | Monoclonal                                                     |
| Isotype               | IgG1                                                           |
| Clone ID              | 1A184                                                          |
| Form                  | Liquid                                                         |
| Storage buffer        | 0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.            |
| Purification          | Affinity chromatography.                                       |
| Immunogen             | Native Human IgM.                                              |
| Stability and Storage | Store at 4°C for 12 months. Protect from light. Do not freeze. |
| Alternative Names     | IgM                                                            |
| Note                  | For research use only.                                         |